FORA 6 Blood Glucose / Hematocrit / Hemoglobin Test Strip

312-0000001-302 ver 1.0 2018/06

Warnings

- For **in vitro** diagnostic use (for use outside of the body only).
- For single use only.
- Healthcare professionals and other users testing multiple patients with this system should handle everything that has come into contact with human blood carefully to prevent transmitting infectious diseases, including sanitized objects.
- Please read this sheet and your FORA 6 Plus / Connect Multi-functional Monitoring System Owner's Manual before you use this test strip. Use only FORA 6 Test Strips with FORA 6 Plus / Connect Multi-functional Monitoring System to obtain accurate results, and be covered by the manufacturer's warranty.
- Results may be inaccurate when testing on patients with abnormally low blood pressure or those who are in shock.
- Please do not use FORA 6 Plus / Connect Multi-functional Monitoring System on critically ill patients. While the blood glucose result is extremely hyperglycemic (over 600 mg/dL (33.3 mmol/L)), the collection of capillary blood from the approved sample sites is not advised when the peripheral circulation is impaired as the delivery of physiological blood glucose level might not be a true reflection. The following circumstances may apply: severe dehydration as a result of diabetic ketoacidosis or due to stress hyperglycemic, hyperosmolar non-ketotic coma, shock, decompensated heart failure NYHA Class IV or peripheral arterial occlusive disease.
- Keep test strips and lancets away from small children. If swallowed, consult a doctor immediately for advice.

Intended Use

FORA 6 test strips, when used together with FORA 6 Plus / Connect Multi-functional Monitoring System is intended for the quantitative measurement of blood glucose, hematocrit and hemoglobin in fresh capillary whole blood from the finger. It is for in vitro diagnostic use only. This system is not intended for use in the diagnosis or screening of diabetes mellitus.

Professionals may test with capillary and venous blood sample; home use is limited to capillary whole blood testing. Use ONLY heparin for anticoagulation of whole blood samples. Please do NOT use EDTA for anticoagulation.

Limitations

- Hematocrit: The hematocrit level is limited to between 0% and 70%. Please ask your healthcare professional if you do not know your hematocrit level.
- Neonatal Use: **This test strip can be used for the testing of newborns.**
- Please see Appendix: Summary of substances and concentrations in excess of limitation with interference.
- Altitude Effects: Altitudes up to 3,275 m (10,742 ft) do not affect test results.
- Do not test blood glucose during or soon after a xylose absorption test. Xylose in the blood can produce elevated glucose results.

Storage and Handling

- Use each test strip immediately after taking it out of the vial or individual foil pack. Close the vial immediately after taking out a strip.
- **IMPORTANT: Do not use the test strips if they have expired.**
- Keep the vial closed at all times.
- Store the test strips in their original vial **ONLY**. Do not transfer them to a new vial or any other containers.
- Do not touch the test strips with wet hands.
- Do not bend, cut, or alter the test strip.
- Store the test strips in a cool, dry place between 2°C and 30°C (35.6°F and 86°F) and below 85% relative humidity.
- Keep the test strips away from direct sunlight. Do not store the test strips in high humidity.

Testing Your Blood Glucose

PLEASE WASH AND DRY YOUR HANDS BEFORE PERFORMING ANY TESTS.

Please refer to your Owner's Manual for more information.

The used lancet and test strip are potentially biohazardous. Please dispose of them carefully according to your local regulations.

Reading Your Result

Your blood glucose readings deliver plasma equivalent results and are displayed in mg/dL or mmol/L. Your hematocrit readings will be displayed in % and hemoglobin readings in g/dL.

Reference values	
Normal plasma glucose range for people without diabetes *1	Fasting and before meal < 100 mg/dL (5.6 mmol/L) 2 hours after meals < 140 mg/dL (7.8 mmol/L)
Hemoglobin *2	Male 14 – 17 g/dL Female 12 – 15 g/dL
Hematocrit *2	Male 41% – 50% Female 36% – 44%

*1. American Diabetes Association (2014). Clinical Practice Recommendations. Diabetes Care, 37 (Supplement 1): S16.
*2. NIH - National Heart, Lung, and Blood Institute (NHLBI) https://www.nhlbi.nih.gov/health-topics/blood-tests

Please consult your doctor to determine a target range that works best for you.

Questionable or inconsistent results
If your test results are unusual or inconsistent with how you are feeling:

- Make sure the confirmation window of the test strip is completely filled with blood.
- Check the expiry date of the test strips.
- Check the performance of your meter and test strip with the control solutions.

Please Note: Unusually high or low blood glucose levels may be symptoms of a serious medical condition. If most of your results are unusually high or low, please contact your healthcare professional.

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Substance	Concentration Tested (mg/dL)	Therapeutic / Physiologic Concentration Range (or Upper Limit) (mg/dL)
Pralidoxime Iodide	> 5,0	~ 10 (i.v. Dose 500 mg)

Referenz / Referenz / Référence / Referenz

- McEnroe, J Robert, et al. National Committee for Clinical Laboratory Standards. Interference testing in clinical chemistry; Approved Guideline - 2nd edition. NCCLS : 2005 - EPT-42, volume 25, number 27.
- EN ISO 15197: 2015 (E): In vitro diagnostic test systems - Requirements for blood-glucose monitoring systems for self-testing in managing diabetes mellitus, Second edition, section 6.4.4 Interference Testing.
- Hardman JG, Limbird LE, Molinoff PB, et al. Goodman & Gilman's the Pharmacological Basis of Therapeutics. 9th ed. New York, NY: McGraw-Hill; 1996: 1570-1571.
- Zaloga GP, Macgregor D. The Critical Care Drug Handbook. 2nd ed. St. Louis, MO: Mosby-Year Book; 1997.
- Kost GJ, Vu HT, Lee JH, et al. Multicenter study of oxygen-intensive handheld glucose point-of-care testing in critical care/hospital/ambulatory patients in the United States and Canada. Crit Care Med. 1998;26:581-590.
- Institut für Diabetes Technologie, System accuracy evaluation of FORA 6 Connect Multi-Functional Monitoring System based on ISO 15197:2013 & EN ISO 15197:2015. Project No.: IDT-1739(2A)-FS; 2018

Appendix / Anhang / Annexe / Appendice

Symbol / Symbol / Symbole / Symbol	Referent / Bedeutung / Signification / Bedeutung
IVD	<ul style="list-style-type: none"> • <i>In vitro</i> diagnostic medical device • Für die <i>In-Vitro</i>-Diagnostik • Dispositif médical de diagnostic <i>in vitro</i> • Dispositivo medico per uso diagnostico <i>in vitro</i>
⊗	<ul style="list-style-type: none"> • Do not reuse • Nicht wiederverwenden • À usage unique • Non riutilizzare
	<ul style="list-style-type: none"> • Consult instructions for use • Gebrauchsanweisung beachten • Consulter le mode d'emploi • Leggere le istruzioni per l'uso
	<ul style="list-style-type: none"> • Temperature limitation • Temperaturbegrenzung • Plage de températures limites • Limite di temperatura
	<ul style="list-style-type: none"> • Use by • Zu verwenden bis • Date limite d'utilisation • Data di scadenza
	<ul style="list-style-type: none"> • Keep Dry • Vor Feuchtigkeit schützen • Gardez à sec • Tenere in luogo asciutto
	<ul style="list-style-type: none"> • Keep away from sunlight • Vor sonnenlicht schützen • Tenir à l'écart de la lumière du soleil • Tenere lontano dai raggi del sole
CE	<ul style="list-style-type: none"> • CE mark • CE-Zeichen • Marquage CE • Marcatura CE
LOT	<ul style="list-style-type: none"> • Batch code • Chargencode • Code du lot • Codice di lotto
	<ul style="list-style-type: none"> • Humidity limitation • Einschränkung bei Luftfeuchtigkeit • Limite d'humidité • Limite di umidità
	<ul style="list-style-type: none"> • Do not use if package is damaged • Bei beschädigter verpackung nicht verwenden • Ne pas utiliser si l'emballage est endommagé • Non utilizzare se la confezione è danneggiata
	<ul style="list-style-type: none"> • Manufacturer • Hersteller • Fabricant • Produttore
REF	<ul style="list-style-type: none"> • Model number • Modell-Nr. • Numéro de modèle • Numero di modello
	<ul style="list-style-type: none"> • Quantity • Menge • Quantité • Quantità

ForaCare Suisse AG
Neugasse 55, 9000,
St. Gallen, Switzerland
www.foracare.ch